



SEALCUFF[®] Waterproof Cast & Bandage Protector

SEALCUFF[®] Waterproof Cast & Bandage Protector is a newest wound care products. It is completely watertight keeping your bandages, prosthesis, and casts dry. Stay dry while taking shower!

Diagrams of products for daily shower or bathing:



Adult Short Arm 330231



Adult Long Arm 330331



Adult Short Leg 330531

Description-Adult (MPE: $\pm 0.5\text{cm}$)

Name	Model	Loop (OD/ID)	Size (LxW)	Box Size (cm)	G.W (kgs)
Adult Hand	330131	15.5/3.5cm	34.5X22cm	19X16X3.5	0.08
Adult Short Arm	330231	15.5/3.5cm	57X23.5cm	19X16X3.5	0.12
Adult Long Arm	330331	19.5/5cm	63X43.5cm	23.5X20X3.5	0.19
Adult Foot/Ankle	330631	19.5/5cm	32.5X35.5cm	23.5X20X3.5	0.14
Adult Short Leg	330531	19.5/5cm	66X43.5cm	23.5X20X3.5	0.20
Adult Long Leg	330431	19.5/5cm	78.5X35.5cm	23.5X20X3.5	0.22
Adult Extra Long Leg	330432	19.5/5cm	104X43.5cm	23.5X20X3.5	0.24
Adult Wide Short Leg	330532	19.5/5cm	61.5X42.5cm	23.5X20X3.5	0.18
Adult Wide Short Arm	330232	19.5/5cm	61X27.5cm	23.5X20X3.5	0.22

Description-Pediatric(Child) (MPE: ±0.5cm)

Name	Model	Loop _(OD/ID)	Size _(LxW)	Box Size _(cm)	G.W _(kgs)
Pediatric Small Arm	330130	15.5/3.5cm	34.5X24cm	19X16X3.5	0.08
Pediatric Medium Arm	330230	15.5/3.5cm	48.5X24cm	19X16X3.5	0.11
Pediatric Large Arm	330330	15.5/3.5cm	53.5X38cm	19X16X3.5	0.12
Pediatric Small Leg	330630	15.5/3.5cm	32X28.5cm	19X16X3.5	0.09
Pediatric Medium Leg	330530	19.5/5cm	44X34.5cm	23.5X20X3.5	0.16
Pediatric Large Leg	330430	19.5/5cm	74.5X44cm	23.5X20X3.5	0.22

Features

- Excellent watertight seal, the material is elastic thermoplastic with high seal for medical use.
- Easy one-handed application, no irritation for skin.
- Non-latex reusable design.
- Adult and pediatric size.

Indication

Ideal for keeping casts, bandages or wounds, IV & PICC lines dry while showering or bathing.

Diagrams of products for daily care



Model:

Size	330734-XS	330733-S	330730-M	330731-L	330732-XL
Circumference	19~22cm	21~24cm	23~26cm	26~29cm	28~31cm

Feature

It is made of high elasticity material, available in breathable, prevent off-line.

Indication

Fixing and protecting the PICC catheter and wound dressing in daily care.